



ENT

SINGLE USE SUCTION TUBE

Product
introduction



CONTENTS

| | | |
|----|--------------------------------------|----|
| 01 | Polished Suctions | 01 |
| 02 | Polished Suctions with Obturator | 02 |
| 03 | Rounded Suctions | 04 |
| 04 | Fenestrated Suctions | 05 |
| 05 | Rounded & Fenestrated Suctions | 06 |
| 06 | Rounded Suction with 4 Fenestrations | 07 |
| 07 | Soft Tip Suction | 08 |
| 08 | Dissector Suction | 08 |



ENT SINGLE USE SUCTION TUBE

ENT Single Use Suction Tube are used in a variety of medical procedures including laryngeal suction, where they irrigate and aspirate the throat, and aural micro-suction, where they are used to remove stubborn or impacted earwax. A range of supporting accessories allow for greater clinical control.


Our core range of ENT Single Use Suction Tube are available in a variety of tube lengths and lumen sizes to suit each case. Selected suction tubes are available with atraumatic tips, with a choice of polished, rounded, fenestrated, rounded tips and so on. Our suction tubes also available for Zoellner, Mastoid, Frazier, Magill.

The Zoellner suction tubes are used in aural micro-suction to remove impacted earwax.

The Mastoid suction tubes are used in radical procedures on the frontal sinuses, in operations on the antrum, in external ethmoidectomy and also in evisceration of the eye and operations on the lacrimal.

The Frazier suction tubes are used to remove fluid or debris from nasal cavities and in laryngeal suction.

The Magill suction tubes are used to remove impacted earwax during aural micro-suction.



Polished Tip

Polished Suctions

All suction tubes have smooth polished tips. Ergonomic balanced design for ease of use. Widely range of lengths and lumens to suit each clinic use.

Purple Handle Is Available



F Features:

- Single use for ENT surgery.
- A variety of tube lengths and lumen sizes for each clinic suit.
- Lightweight, non-blocking design.
- Ergonomic design polymer handle comfortably suits left or right-handed users.
- Available with control hole handle, allows greater suction when covered with user's finger.

| Product Code | Size | Description |
|--------------|---------------------------------|-----------------------------|
| 1262S20W1 | Working Length 85/60mm, OD=2mm | With Control vent, Zoellner |
| 1262S30W1 | Working Length 85/60mm, OD=3mm | With Control vent, Mastoid |
| 1262S20W2 | Working Length 128/90mm, OD=2mm | With Control vent, Frazier |
| 1262S30W2 | Working Length 128/90mm, OD=3mm | With Control vent, Frazier |


| | | |
|-----------|----------------------------------|--|
| 1262S40W2 | Working Length 128/90mm, OD=4mm | With Control vent, Frazier |
| 1262S50W2 | Working Length 128/90mm, OD=5mm | With Control vent, Frazier |
| 1262S40W3 | Working Length 138/100mm, OD=4mm | With Control vent, Frazier |
| 1262S20W4 | Working Length 168/130mm, OD=2mm | With Control vent, Frazier |
| 1262S30W4 | Working Length 168/130mm, OD=3mm | With Control vent, Frazier |
| 1262S40W4 | Working Length 168/130mm, OD=4mm | With Control vent, Frazier |
| 1262S50W4 | Working Length 168/130mm, OD=5mm | With Control vent, Frazier |
| 1262R20W0 | Working Length 85/35mm, OD=2mm | With Control vent, Magill |
| 1262R30W0 | Working Length 85/35mm, OD=3mm | With Control vent, Magill |
| 1264R20W0 | Working Length 85/35mm, OD=2mm | Purple handle With Polished Tip, Frazier |
| 1264S20W1 | Working Length 85/60mm, OD=2mm | Purple handle With Polished Tip, Frazier |
| 1264S30W2 | Working Length 128/90mm, OD=3mm | Purple handle With Polished Tip, Frazier |



Polished Suctions with Obturator

F Features:

- Single use for ENT surgery.
- Flexible cannula could be bent to needed position without occlusion.
- Available with integrated transparent handle for visual check.



Polished Tip with obturator

| Product Code | Size | Description |
|--------------|------|--|
| 1242A08 | Fr08 | With flexible cannula and obturator, Frazier |
| 1242A10 | Fr10 | With flexible cannula and obturator, Frazier |
| 1242A12 | Fr12 | With flexible cannula and obturator, Frazier |
| 1242A18 | Fr18 | With flexible cannula and obturator, Frazier |

Transparent Handle Polished Suctions With Obturator



| Product Code | Size | Description |
|--------------|------|---|
| 1241A06 | Fr06 | Transparent handle, With flexible cannula and obturator, Frazier |
| 1241A08 | Fr08 | Transparent handle, With flexible cannula and obturator, Frazier |
| 1241A10 | Fr10 | Transparent handle, With flexible cannula and obturator, Frazier |
| 1241A12 | Fr12 | Transparent handle, With flexible cannula and obturator, Frazier |
| 1241B06 | Fr06 | Transparent handle, With flexible cannula and X-ray line obturator, Frazier |
| 1241B08 | Fr08 | Transparent handle, With flexible cannula and X-ray line obturator, Frazier |
| 1241B10 | Fr10 | Transparent handle, With flexible cannula and X-ray line obturator, Frazier |
| 1241B12 | Fr12 | Transparent handle, With flexible cannula and X-ray line obturator, Frazier |



Rounded Suctions

A range of suction tubes with a smooth rounded tip for atraumatic tissue contact. Used for ENT (Ear-Nose-Throat) surgery.



F Features:

- Fine smooth rounded tip for atraumatic tissue contact.
- Available with a variety of lengths and lumens for each cases.




Rounded Tip

| Product Code | Size | Description |
|--------------|----------------------------------|---------------------------|
| 1262P40W1 | Working Length 85/60mm, OD=4mm | With Rounded Tip, Frazier |
| 1262P20W2 | Working Length 128/90mm, OD=2mm | With Rounded Tip, Frazier |
| 1262P30W2 | Working Length 128/90mm, OD=3mm | With Rounded Tip, Frazier |
| 1262P40W2 | Working Length 128/90mm, OD=4mm | With Rounded Tip, Frazier |
| 1262P20W4 | Working Length 168/130mm, OD=2mm | With Rounded Tip, Frazier |
| 1262P30W4 | Working Length 168/130mm, OD=3mm | With Rounded Tip, Frazier |
| 1262P40W4 | Working Length 168/130mm, OD=4mm | With Rounded Tip, Frazier |
| 1262O30W0 | Working Length 85/35mm, OD=3mm | With Rounded Tip, Magill |



Fenestrated Suctions

A range of suction with fenestrated tips to disperse aspiration and prevent suction trauma.



Fenestrated Tip

| Product Code | Size | Description |
|--------------|----------------------------------|-------------------------------|
| 1262T20W1 | Working Length 85/60mm, OD=2mm | With Fenestrated Tip, Frazier |
| 1262T30W1 | Working Length 85/60mm, OD=3mm | With Fenestrated Tip, Frazier |
| 1262T40W1 | Working Length 85/60mm, OD=4mm | With Fenestrated Tip, Frazier |
| 1262T20W2 | Working Length 128/90mm, OD=2mm | With Fenestrated Tip, Frazier |
| 1262T30W2 | Working Length 128/90mm, OD=3mm | With Fenestrated Tip, Frazier |
| 1262T20W4 | Working Length 168/130mm, OD=2mm | With Fenestrated Tip, Frazier |
| 1262T30W4 | Working Length 168/130mm, OD=3mm | With Fenestrated Tip, Frazier |
| 1262T40W4 | Working Length 168/130mm, OD=4mm | With Fenestrated Tip, Frazier |



Rounded & Fenestrated Suctions

Suctions that benefit from a tip that is both fenestrated and rounded.

F Features:

- Suction control.
- For micro-laryngoscopy (MLS).
- Rounded opening for suction protection.



Rounded & Fenestrated Tip

| Product Code | Size | Description |
|--------------|----------------------------------|---|
| 1262Q20W2 | Working Length 128/90mm, OD=2mm | With Rounded & Fenestrated Tip, Frazier |
| 1262Q30W2 | Working Length 128/90mm, OD=3mm | With Rounded & Fenestrated Tip, Frazier |
| 1262Q40W2 | Working Length 128/90mm, OD=4mm | With Rounded & Fenestrated Tip, Frazier |
| 1262Q20W4 | Working Length 168/130mm, OD=2mm | With Rounded & Fenestrated Tip, Frazier |
| 1262Q30W4 | Working Length 168/130mm, OD=3mm | With Rounded & Fenestrated Tip, Frazier |
| 1262Q40W4 | Working Length 168/130mm, OD=4mm | With Rounded & Fenestrated Tip, Frazier |




Fenestrated Suctions

This ultra-delicate suction has 4 fenestrations that are 1mm in diameter at the tip.



Rounded Suction with 4 Fenestrations Tip

| Product Code | Size | Description |
|--------------|---------------------------------|--|
| 1262H40W2 | Working Length 128/90mm, OD=4mm | Rounded Tip, Fenestrations 4*1 mm, Mastoid |



Soft Tip



Soft Tip Suction


A suction with a soft silicone tip for very delicate procedures.

| Product Code | Size | Description |
|--------------|---------------------------------|---------------|
| 1262M20W2 | Working Length 128/90mm, OD=2mm | With Soft Tip |



Dissector Suction

A straight tube finished with a flattened polished tip.



Dissector Flattened Tip



| Product Code | Size | Description |
|--------------|----------------------------------|------------------------------|
| 1262N40W4 | Working Length 168/130mm, OD=4mm | With Dissector Flattened Tip |



ENT

SINGLE USE SUCTION
TUBE

CONTACT US



0086-510-88880366



LOUJIN INDUSTRIAL PARK, SHUOFANG,
WUXI, 214143, P.R.CHINA



INFO@GCMEDICA.COM



WWW.GCMEDICA.COM
